

IDENTIFYING LEAKS IN PLASTIC ANTIFREEZE CONTAINERS

Tested: Plastic Antifreeze Bottles

Tested with: TapTone 1000-C Compression System

The purpose of this test was to find leaks in plastic antifreeze bottles using the TapTone 1000-C compression system. Leaking antifreeze bottles can be hazardous for consumers since they contain harsh chemicals.

The T1000-C Compression Sensor is ideal for finding potential leakers in antifreeze containers before they leave the processing plant.



Plastic antifreeze bottle



Bottle being tested with the T1000-C

How It Works

The T1000-C Compression Sensor utilizes dual parallel belts to apply force to the sidewall of a passing container. This action also compresses the headspace of the container which allows a load cell to take a force measurement at the discharge of the system. Utilizing DSP technology, the controller analyzes the measurement and assigns a merit value to each container. If the merit value is outside of the acceptable range, a reject signal activates a remote reject system.

See next page for test results



T-1000 Compression Controller and Sensor. Sensor has a cantilever design that suspends over the existing conveyor.



TELEDYNE TAPTONE
Package Inspection Systems
A Teledyne Technologies Company


www.TapTone.com

TEST

Empty flexible plastic containers were supplied for this evaluation. The plastic twist-off caps have an induction seal lining that is affixed to the containers after filling. For testing purposes, leaks of three different sizes were evaluated: 0.51 mm, 0.64 mm, 0.76 mm (0.020, 0.025 and 0.030 inches). The plastic bottles were oriented to pass through the compression belts with the sides of the container against the belt, allowing for a wider surface area against the compression sensor.



During the compression cycle, the sealed, non-leaking containers generated an average merit value* of 875. When leaking containers were introduced, those with hole sizes of 0.51 mm, 0.64 mm and 0.76 mm (0.020, 0.025 and 0.030 inches) produced average merit values of 500, 425, and 300 respectively.



*Merit value is a calculated number determined using an algorithm to compute a resultant from a set of data values.

Note: It is possible for the plastic cap to be torqued tight enough to create a seal on the container even if the foil seal inside is leaking. In this case the system would consider the container to be a non-leaker.

SUMMARY

Testing was successful with the T1000-C. In the samples provided, the system was capable of consistently detecting leaks with a diameter as small as 0.51mm (0.020 inch).



TELEDYNE TAPTONE
Package Inspection Systems
 A Teledyne Technologies Company

www.TapTone.com

49 Edgerton Drive • North Falmouth, MA 02556 USA

P: 508.563.1000

F: 508.564.9945

E-Mail: taptone@teledyne.com

ISO 9001 Certified CE Compliant

Specifications subject to change without notice.
 TAPTONE is a registered trademark of TELEDYNE TAPTONE. Copyright 2007, TELEDYNE TAPTONE.